



Nanosorter: A high-throughput nanoshell sorting device

Team Tiny Tools, Dept. of Bioengineering, Rice University
 Drew Bryant, Dania Daye, Joshua Carroll, and Adrielle Fry
 {dbryant,dania,joshc,afry}@rice.edu



Motivation


Nanoshells offer promising oncology treatment solutions

- Non-invasive destruction of cancer tissue
- Low collateral damage to surrounding healthy tissue
- Superior biocompatibility and side effect profile
- Treatment individualized to patients

Problem:

- Nanoshell production methods yield unwanted by-products
- Current purification methods cause a manufacturing bottleneck
- Current by-product removal techniques are low yield and time-consuming
- Size exclusion chromatography recovers less than 1% of input nanoshells
- Repeated centrifugation recovers approximately 30% of input nanoshells

Crude Nanoshell Solution



Solution:

- Nanosorter implements a high-throughput, scalable nanoshell purification process
- Generalizable to other types of nanoparticles


Design Objectives

- High-throughput, scalable process
 - Repeatable nanoshell purification
 - Greater than 40% recovery of desired nanoshells input (> 40% yield)
 - Per device production cost < \$2000
 - Device footprint < 4 ft²

Acknowledgements


We would like to thank André Gobin, Dr. Maria Oden, Dr. Beth Beason, The Brown Foundation Teaching Grant, the Rice Bioengineering Dept., and funding from the Center for Biological and Environmental Nanotechnology.

Protein-to-nanoshell Conjugation




- Bifunctional poly-ethylene glycol (Bi-PEG) links proteins to nanoshell surface
- Remaining nanoshell surface area coated with PEG-thiol (PEG-SH)
 - Reduces nanoshell aggregation
 - Promotes longevity of nanoshells

Size-Selective Nanoshell Heating



- Nanoshell light absorption wavelength is governed by nanoshell size and shell thickness
- Applying high intensity laser to a sample causes some nanoshells to heat up
- Nanoshell heat radiation denatures (unfolds) conjugated protein
- Only nanoshells of resonance corresponding to the applied laser wavelength become hot

Hydrophobic Interaction Chromatography



- Desired nanoparticles have denatured (unfolded) protein attached
 - Denatured proteins have many exposed hydrophobic regions
 - Hydrophobic regions have strong attraction to hydrophobic beads
- Undesired particles have native (folded) protein attached
 - Very few exposed hydrophobic regions
 - Weak attraction to hydrophobic beads

The Nanosorter Device

Input:


Crude nanoshell solution where both nanoshells and by-products are protein conjugated

Purification:

Input nanoshell solution circulated through heating stage multiple times
 Heat-treated nanoshell solution sent through hydrophobic interaction column

Output:


Multiple solutions of fractionated nanoshell product



Result: removal of unwanted by-product

- Spectral analysis of crude nanoshell sample before nanosorting shows large amount of undesired by-product.
- After passing crude sample through nanosorter, a relative decrease in by-product population is shown.
- Nanosorter does not produce 100% yield of desired nanoshells, but ability to purify crude samples is demonstrated.

- Successful proof of concept



Conclusions

The size-selective purification method presented has shown to be a viable alternative to current nanoshell separation techniques.

Future Work

Examine the effect of adjusting experimental parameters to optimally enhance hydrophobic interactions during chromatography stage

Optimize size-selective nanoshell heating running time to maximally denature protein conjugated to desired nanoshells while minimizing non-specific heating